



Mouse Anti-19-methyltestosterone monoclonal antibody, clone NT (CABT-L3097)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Methyltestosterone with carrier protein.
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	NT
Purification	Protein A purified.
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice
Warnings	For research use only.

BACKGROUND

Introduction Methyltestosterone and methylboldenone are orally active synthetic derivatives of testosterone. In addition to human clinical applications, and illegal use as a growth promoter in cattle, methyltestosterone has been used as a masculinising agent in aquaculture. Methyltestosterone and methylboldenone are mainly excreted in urine as glucuronides of hydroxylated metabolites.

Keywords

Methyltestosterone;Android;Testred;Virilon;17b-hydroxy-17a-methylandro-4-ene-3-one;METANDREN;MESTERONE;3H-Phenothiazine, 7-amino-3-imino, monohydrochloride;Methyl-1-testosterone;17ALPHA-METHYLTESTOSTERONE
