



Methyl Mercury Chloride [BSA] (DAG-P2514B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Nature	Synthetic
Expression System	N/A
Conjugate	BSA
Applications	ELISA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1mg
Buffer	0.01M pH7.4 PBS
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.