



Clopidol [BSA] (DAG010S)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	N/A
Nature	Synthetic
Expression System	N/A
Species	N/A
Conjugate	BSA
Applications	ELISA, LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M PBS, pH7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Ship	Wet ice
Warnings	This product is for research use only and is not intended for diagnostic use.

BACKGROUND

Introduction Clopidol is an organic compound that is used as in veterinary medicine as a coccidiostat. It is prepared industrially by a multistep process from dehydroacetic acid.

Keywords Clopidol