



# Clopidol [KLH] (DAG012S)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	N/A
<b>Nature</b>	Synthetic
<b>Expression System</b>	N/A
<b>Species</b>	N/A
<b>Conjugate</b>	KLH
<b>Applications</b>	ELISA, LFIA
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1 mg
<b>Buffer</b>	0.01M PBS, pH7.4
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
<b>Ship</b>	Wet ice
<b>Warnings</b>	This product is for research use only and is not intended for diagnostic use.

## BACKGROUND

**Introduction** Clopidol is an organic compound that is used as in veterinary medicine as a coccidiostat. It is prepared industrially by a multistep process from dehydroacetic acid.

**Keywords** Clopidol