



Tobramycin [BSA] (DAG3141)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Tobramycin, BSA-Conjugated
Specificity	14/1 (Moles Antigen/Moles Protein)
Nature	Synthetic
Expression System	N/A
Species	N/A
Conjugate	BSA
Applications	immunogen
Procedure	None
Preservative	None
Storage	2-8°C short term, -20°C long term
Warnings	PLEASE note that this product is intended for research use only; not for diagnostic or clinical use.

BACKGROUND

Introduction	Tobramycin is an aminoglycoside antibiotic used to treat infections caused by many different bacteria. Tobramycin is usually administered by intravenous infusion, intramuscular injection, or inhalation.
Keywords	Tobramycin