



## Robenidine [OVA] (DAG5621)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

**Antigen Description** Robenidine is a coccidiostat. It is an antibiotic used for the control of coccidiosis, a debilitating protozoal infection in poultry. Although there are alternative antibiotics available, robenidine is important in the management of antibiotic resistance as farmers rotate the use of robenidine with other antibiotics to try to preserve the effectiveness of these products in fighting infections.

<b>Nature</b>	Synthetic
<b>Expression System</b>	N/A
<b>Species</b>	N/A
<b>Conjugate</b>	OVA
<b>Applications</b>	The antigen can be used in ELISA assay.
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1 mg
<b>Buffer</b>	0.01M pH7.4 PBS
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.