



Pb2+ [KLH] (DAGA-005K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Nature	Synthetic
Expression System	N/A
Conjugate	KLH
Applications	Immunogen
Procedure	None
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords Lead(2+);Lead (II) ion;Lead(2+) ion;Lead ion;Omaha;Pb2+