



Cd²⁺ [HRP] (DAGA-006H)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Nature	Synthetic
Expression System	N/A
Conjugate	HRP
Applications	ELISA
Procedure	None
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords Cadmium ion;Cadmium(2+);Cadmium(2+) ion;Cadmium cation;Cadmium ions;Cadmium Atom