



Recombinant DENV type 1 Non-structural Protein 2B (Fragment) [His] (DAGA-1045)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|----------------------------|--|
| Antigen Description | Recombinant 6xHis tagged DENV type 1 Non-structural protein 2B (Fragment) (UniProtKB No. A0A1L2BPD7) was expressed. |
| Nature | Recombinant |
| Expression System | E. coli |
| Purity | >85% |
| Conjugate | His |
| Applications | ELISA, Immunogen |
| Recommended Usage | This item requires custom production and lead time is between 4-8 weeks. We can custom produce based on your specifications. |
| Procedure | None |
| Format | Lyophilized |
| Concentration | Batch dependent - please inquire should you have specific requirements. |
| Buffer | Each vial contains 100 µg of purified protein in PBS |
| Preservative | None |
| Storage | Store at 2-8 °C short term and -20°C long term. It is recommended to aliquot needed volumes and store at -20°C. Avoid repeated freeze/thaw cycles. |

BACKGROUND

Keywords DENV type 1;Dengue Virus type 1;Dengue Virus 1+2+3+4;Dengue Virus Type 1, 2, 3, 4;Group IV;Flaviviridae;Flavivirus;Dengue virus;Dengue Virus;DENV

ANTIGEN GENE INFORMATION

