



## Acrylamide [KLH] (DAGA-110K)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Nature</b>	Synthetic
<b>Expression System</b>	N/A
<b>Conjugate</b>	KLH
<b>Applications</b>	Immunogen
<b>Procedure</b>	None
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1mg
<b>Buffer</b>	0.01M pH7.4 PBS
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### BACKGROUND

<b>Keywords</b>	Acrylamide
-----------------	------------