



Vitamin B5 [HRP] (DAGA-123H)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Pantothenic acid; pantothenate; vitamin B5; Pantothenic; Ca pantothenaat; kyselínápantothenová; d-Pantothenic Acid; Vitamin B5
Nature	Synthetic
Expression System	N/A
Conjugate	HRP
Applications	ELISA
Procedure	None
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	D-pantothenic acid; Pantothenic acid; Vitamin B5;
-----------------	---