



# Vitamin B5 [KLH] (DAGA-123K)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Antigen Description</b>	Pantothenic acid; pantothenate; vitamin B5; Pantothenic; Ca pantothenaat; kyselina pantothenova; d-Pantothenic Acid; Vitamin B5
<b>Nature</b>	Synthetic
<b>Expression System</b>	N/A
<b>Conjugate</b>	KLH
<b>Applications</b>	Immunogen
<b>Procedure</b>	None
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1mg
<b>Buffer</b>	0.01M pH7.4 PBS
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## BACKGROUND

<b>Keywords</b>	D-pantothenic acid; Pantothenic acid; Vitamin B5;
-----------------	---