



Recombinant Chikungunya E1 E2 Composite Antigen (DAGA-163)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description Chikungunya is an infection caused by the chikungunya virus. The disease features the sudden onset of fever two to four days after exposure. The fever usually lasts two to seven days, while accompanying joint pains typically last weeks or months but sometimes years. The mortality rate is a little less than 1 in 1,000; the elderly or those with underlying chronic medical problems are most likely to have severe complications.

Nature	Recombinant
Expression System	E. coli
Conjugate	Unconjugated
Applications	WB standard, antibody ELISA, immunogen, etc.
Procedure	None
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1mg, 10mg
Preservative	None

BACKGROUND

Keywords CHIKV;Chikungunya Virus;Chikungunya E1;Chikungunya E2;CHIKV E1;CHIKV E2