



# Trenbolone [BSA] (DAGA-178B)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

**Antigen Description** Trenbolone is a steroid used by veterinarians on livestock to increase muscle growth and appetite. To increase its effective half-life, trenbolone is not used in an unrefined form, but is rather administered as ester derivatives such as trenbolone acetate

**Nature** Synthetic

**Expression System** N/A

**Conjugate** BSA

**Applications** LFIA

**Procedure** None

**Format** Liquid

**Concentration** Batch dependent - please inquire should you have specific requirements.

**Size** 1mg

**Buffer** 0.01M pH7.4 PBS

**Preservative** None

**Storage** Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## BACKGROUND

**Keywords** Trenbolone; 17 $\beta$ -Hydroxyestra-4; 9; 11-trien-3-one