



## Tobramycin [KLH] (DAGA-213K)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Antigen Description</b>	Tobramycin is an aminoglycoside antibiotic used to treat various types of bacterial infections, particularly Gram-negative infections.
<b>Nature</b>	Synthetic
<b>Expression System</b>	N/A
<b>Conjugate</b>	KLH
<b>Applications</b>	Immunogen
<b>Procedure</b>	None
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1mg
<b>Buffer</b>	0.01M pH7.4 PBS
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### BACKGROUND

**Keywords** Tobramycin; Anti – Tobramycin; 1-epitobramycin; 3'-deoxykanamycinb; d-6-tyrideoxy-alpha-d-ribohexopyranosyl-(1-6))-2-deoxy; distobram gernebcin lpha-d-ribohexopyranosyl-(1-4))-2-deoxy-d-streptamine; nebcin; nebicina; nebramycin6; nebramycinvi