



## Decoquinatate [BSA] (DAGA-249B)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Antigen Description</b>	Decoquinatate is a quinolone anticoccidial which is nontoxic but subject to a high rate of development of resistant coccidia.
<b>Nature</b>	Synthetic
<b>Expression System</b>	N/A
<b>Conjugate</b>	BSA
<b>Applications</b>	LFIA
<b>Procedure</b>	None
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1mg
<b>Buffer</b>	0.01M pH7.4 PBS
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### BACKGROUND

<b>Keywords</b>	Decoquinatate; Coccidiostat;
-----------------	------------------------------