



Decoquinatate [KLH] (DAGA-249K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Decoquinatate is a quinolone anticoccidial which is nontoxic but subject to a high rate of development of resistant coccidia.
Nature	Synthetic
Expression System	N/A
Conjugate	KLH
Applications	Immunogen
Procedure	None
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	Decoquinatate; Coccidiostat;
-----------------	------------------------------