



Lasalocid [KLH] (DAGA-282K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Lasalocid is able to make neutral complexes with monovalent and divalent cations and transport them through apolar phase (including lipid bilayer membranes). It can transport also big organic cations like dopamine.
Nature	Synthetic
Expression System	N/A
Conjugate	KLH
Applications	Immunogen
Procedure	None
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords Lasalocid; Bovatec;