



Maduramicin [KLH] (DAGA-290K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Maduramicin is a monoglycoside polyether derived from the fungus <i>Actinomadura rubra</i> . It is classified as an ionophor and predominantly used as an anti-coccidial agent, which is mixed with poultry feed in the ratio of 6ppm.
Nature	Synthetic
Expression System	N/A
Conjugate	KLH
Applications	Immunogen
Procedure	None
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords Maduramicin; Maduramicin A