



CHIKV Capsid protein [His] (DAGA-292)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	A member of the Togaviridae family, and Alphavirus genus, belonging to the Semliki Forest serological complex. CHIKV is a spherical enveloped virion that measures 60 to 70 nm in diameter, and contains a single-stranded, positive-sense RNA genome.
Nature	Recombinant
Expression System	E. coli
Purity	>90% , based on SDS PAGE
Conjugate	His
Applications	WB standard, antibody ELISA, immunogen, etc.
Procedure	None
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	100 µg, 1 mg
Buffer	In PBS with 8M Urea
Preservative	None
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.

BACKGROUND

Keywords Chikungunya virus capsid protein;CHIKV capsid protein;CHIKV C protein;CHIKV