



Recombinant C. Jejuni Omp protein [His] (DAGA-3031)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Campylobacter jejuni Outer Membrane Protein (OMP) Recombinant Product contains fusion partner His tag.
Nature	Recombinant
Expression System	E. coli
Purity	Multi-step procedure including affinity chromatography.
Conjugate	His
Applications	Suitable for use as control in Lateral Flow Assays. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Molecular Weight	18 kDa
Format	Purified
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	0.1 MG
Buffer	Tris Buffer
Preservative	0.09% Sodium Azide
Storage	Store at 2-8°C for < 1 week. Long term storage at -20°C. Aliquot to avoid multiple freeze/thaw cycles.