



Recombinant E. coli O157 VT1 [His] (DAGA-3082)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Nature	Recombinant
Expression System	E. coli
Purity	Multi-step procedure including affinity chromatography.
Conjugate	His
Applications	Suitable for use as a control in Lateral Flow Assays and as an immunogen. Suitable for spiking into stool samples. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	0.1 MG
Buffer	Carbonate Buffer with salts, pH 8.3.
Preservative	0.09% Sodium Azide
Storage	Store at 2-8°C.