



Recombinant E. coli O157 VT2 [His] (DAGA-3083)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|--------------------------|--|
| Nature | Recombinant |
| Expression System | E. coli |
| Purity | Multi-step procedure including affinity chromatography. |
| Conjugate | His |
| Applications | Suitable for use as a control in Lateral Flow Assays and as an immunogen. Suitable for spiking into stool samples. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded. |
| Format | Purified |
| Concentration | Batch dependent - please inquire should you have specific requirements. |
| Size | 0.1 MG |
| Buffer | Carbonate Buffer with salts, pH 8.3. |
| Preservative | 0.09% Sodium Azide |
| Storage | Store at 2-8°C. |