



Robenidine [KLH] (DAGA-328K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Robenidine is an antibiotic used for the control of coccidiosis, a debilitating protozoal infection in poultry.
Nature	Synthetic
Expression System	N/A
Conjugate	KLH
Applications	Immunogen
Procedure	None
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords Robenidine; Robenidine hydrochloride; Robenzidene;