



Human SAA reference serum (DAGA-652)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Serum amyloid A (SAA) proteins are a family of apolipoproteins associated with high-density lipoprotein (HDL) in plasma. Different isoforms of SAA are expressed constitutively (constitutive SAAs) at different levels or in response to inflammatory stimuli (acute phase SAAs). These proteins are produced predominantly by the liver. The conservation of these proteins throughout invertebrates and vertebrates suggests that SAAs play a highly essential role in all animals.
Specificity	Human
Nature	Native
Expression System	Human
Conjugate	Unconjugated
Procedure	None
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1ml
Preservative	0.1% Sodium Azide
Storage	Frozen -20°C

BACKGROUND

Keywords Human;SAA;Serum

ANTIGEN GENE INFORMATION

MIM None