



Rat cystatin C reference serum (DAGA-689)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description Cystatin C or cystatin 3 (formerly gamma trace, post-gamma-globulin, or neuroendocrine basic polypeptide), a protein encoded by the CST3 gene, is mainly used as a biomarker of kidney function. It has been studied for its role in predicting new-onset or deteriorating cardiovascular disease. It also seems to play a role in brain disorders involving amyloid (a specific type of protein deposition), such as Alzheimer's disease.

Specificity	Rat
Nature	Native
Expression System	Rat
Conjugate	Unconjugated
Procedure	None
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1ml
Preservative	0.1% Sodium Azide
Storage	Frozen -20°C

BACKGROUND

Keywords Rat;Cystatin C;Serum

ANTIGEN GENE INFORMATION

MIM None