



Rat hemopexin reference serum (DAGA-694)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description Hemopexin (or haemopexin; HPX), also known as beta-1B-glycoprotein is a protein that in humans is encoded by the HPX gene and belongs to hemopexin family of proteins. Hemoglobin and its scavenger protein hemopexin (Hx) associate with HDL and influence the inflammatory properties of HDL. In addition it can also be said that HDL from Hx-null mice is proinflammatory. Moreover, Hemopexin deficiency is associated with various other inflammatory diseases such as septic shock and experimental autoimmune encephalomyelitis.

Specificity	Rat
Nature	Native
Expression System	Rat
Conjugate	Unconjugated
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1ml
Preservative	0.1% sodium azide
Storage	Frozen -20°C

BACKGROUND

Keywords Rat;Hemopexin;Serum