



Rat L-FABP reference serum (DAGA-696)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description FABP1 is a human gene coding for the protein product FABP1 (Fatty Acid-Binding Protein 1). It is also frequently known as liver-type fatty acid-binding protein (LFABP). FABP1 is primarily expressed in the liver where it is involved in the binding, transport and metabolism of long-chain fatty acids (LCFAs), endocannabinoids, and other hydrophobic molecules.[3][4][5] Altered expression of the protein has been linked to metabolic conditions including obesity.

Specificity	Rat
Nature	Native
Expression System	Rat
Conjugate	Unconjugated
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1ml
Preservative	0.1% sodium azide
Storage	Frozen -20°C

BACKGROUND

Keywords Rat;L-FABP;Serum