



## Rat prealbumin reference serum (DAGA-700)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Antigen Description</b>	Prealbumin (PAB) is a maker of nutritional status and inflammation of potential prognostic value in acute heart failure (HF).PAB can be used as an indicator of liver damage caused by asphyxia in neonates, and can be used to assess the degree of asphyxia
<b>Specificity</b>	Rat
<b>Nature</b>	Native
<b>Expression System</b>	Rat
<b>Conjugate</b>	Unconjugated
<b>Procedure</b>	None
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements
<b>Size</b>	1ml
<b>Preservative</b>	0.1% Sodium Azide
<b>Storage</b>	Frozen -20°C

### BACKGROUND

**Keywords** Rat;Prealbumin;Serum