



# Rat retinol binding protein reference serum (DAGA-701)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Antigen Description</b>	Retinol-binding proteins (RBP) are a family of proteins with diverse functions. They are carrier proteins that bind retinol. Assessment of retinol-binding protein is used to determine visceral protein mass in health-related nutritional studies.
<b>Specificity</b>	Rat
<b>Nature</b>	Native
<b>Expression System</b>	Rat
<b>Conjugate</b>	Unconjugated
<b>Procedure</b>	None
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements
<b>Size</b>	1ml
<b>Preservative</b>	0.1% Sodium Azide
<b>Storage</b>	Frozen -20°C

## BACKGROUND

**Keywords** Rat;Retinol Binding Protein;Serum