



# Rat vitamin D binding protein reference serum (DAGA-705)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

**Antigen Description** Vitamin D-binding protein, also/originally known as gc-globulin (group-specific component), is a protein that in humans is encoded by the GC gene. Vitamin D-binding protein belongs to the albumin gene family, together with human serum albumin and alpha-fetoprotein. It is a multifunctional protein found in plasma, ascitic fluid, cerebrospinal fluid and on the surface of many cell types. It binds to vitamin D and its plasma metabolites and transports them to target tissues.

<b>Specificity</b>	Rat
<b>Nature</b>	Native
<b>Expression System</b>	Rat
<b>Conjugate</b>	Unconjugated
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements
<b>Size</b>	1ml
<b>Preservative</b>	0.1% sodium azide
<b>Storage</b>	Frozen -20°C

## BACKGROUND

**Keywords** Rat;Vitamin D Binding Protein;Serum