



V. cholerae Cholera Toxin B [FITC] (DAGB148)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Nature	Native
Expression System	V. cholerae
Species	V. cholerae
Conjugate	FITC
Applications	ID, HA
Reconstitution	When reconstituted with 0.4 ml of water, each vial contains 0.2 mg of protein in 0.01M sodium phosphate, pH 7.5.
Procedure	None
Format	Lyophilized power, sealed under vacuum
Size	0.2 mg
Preservative	None
Storage	Store at 4°C prior to and following reconstitution. This product is light sensitive.