



## V. cholerae Cholera Toxin B [Biotin] (DAGB151)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Nature</b>	Native
<b>Expression System</b>	V. cholerae
<b>Species</b>	V. cholerae
<b>Conjugate</b>	Biotin
<b>Applications</b>	ELISA
<b>Reconstitution</b>	When reconstituted with 0.4 ml of water the protein is in 0.01M sodium phosphate, pH7.5.
<b>Endotoxin</b>	The endotoxin content, determined using a kinetic chromagenic LAL assay, is <10 EU/mg.
<b>Procedure</b>	None
<b>Format</b>	Lyophilized power, sealed under vacuum
<b>Size</b>	0.2 mg
<b>Preservative</b>	None
<b>Storage</b>	Store at 2-8°C prior to and after reconstitution.