



Cynomolgus Monkey Sterile Filtered Serum (DAGC284)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview The material for this product was processed in the USA from the blood of healthy Cynomolgus monkeys that had not been inoculated with or exposed to any infectious agents. At no time during collection or processing was the material commingled with any other material of animal origin.

This product is not to be used for direct or indirect exposure or inoculation into laboratory and domestic animals.

Nature	Native
Expression System	Cynomolgus Monkey serum, mixed sex, mixed breed
Species	Cynomolgus Monkey
Purity	0.2 micron (Absolute)
Conjugate	Unconjugated
Applications	ELISA
Procedure	None
Format	Liquid
Concentration	Hemoglobin: <30 mg/dl Protein:>4.0 g/dl
Size	50ml, 100ml, 500ml
Buffer	pH: 6.8 - 8.2
Preservative	None
Storage	This serum is stable for ≥ 3 years when stored -20°C ± 5°C. Avoid repeated freezing and thawing.

BACKGROUND

Keywords

Cynomolgus Monkey Serum; Serum; Monkey Serum
