



Testosterone-3 Antigen [BSA] (DAGF-024)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description Testosterone is the primary male sex hormone and an anabolic steroid. In male humans, testosterone plays a key role in the development of male reproductive tissues such as the testis and prostate, as well as promoting secondary sexual characteristics such as increased muscle and bone mass, and the growth of body hair.

Nature	Synthetic
Expression System	Synthetic
Conjugate	BSA
Procedure	None
Format	Liquid
Size	1mg
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Preservative	0.01% Sodium Azide
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.