



Recombinant HPV16 E6 protein [His] (DAGF-094)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Recombinant Papillomavirus 16 E6 antigen produced in E.Coli contains 159 amino covering the full length of HPV-16 having its C-terminus fused with a 6xHis tag and a total molecular weight of 21kDa. The HPV16 was purified by proprietary chromatographic technique.
Nature	Recombinant
Expression System	E. coli
Purity	>95%
Conjugate	His
Molecular Weight	21kDa
Procedure	None
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Buffer	PBS and 4M Urea
Preservative	None
Storage	Recombinant HPV-16 E6 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.

BACKGROUND

Keywords HPV 16; HPV; Human papillomavirus (HPV); Papillomaviridae; Human Papilloma Virus Type 16; Human Papilloma Virus; Human papillomavirus type 16 E6; Human papillomavirus type 16 E6

ANTIGEN GENE INFORMATION

MIM None