



Nicarbazine standard (DAGS084)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Use this reference material as calibrant for chromatography or any other analytical technique.
Expression System	N/A
Purity	98.79%(Please refer to the certificate of analysis of each batch.)
Conjugate	N/A
Cas.No	330-95-0
Molecular Weight	426.38g/mol
Format	Neat
Size	100 mg
Storage	The reference material should be stored in the original sealed container at the 20±4°C. After use the bottle should be tightly closed and protected from moisture.
Ship	Ship at room temperature.