



Anti-CAMKK2 monoclonal antibody (DCABH-10825)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description The product of this gene belongs to the Serine/Threonine protein kinase family, and to the Ca(2+)/calmodulin-dependent protein kinase subfamily. This protein plays a role in the calcium/calmodulin-dependent (CaM) kinase cascade by phosphorylating the downstream kinases CaMK1 and CaMK4. Seven transcript variants encoding six distinct isoforms have been identified for this gene. Additional splice variants have been described but their full-length nature has not been determined. The identified isoforms exhibit a distinct ability to undergo autophosphorylation and to phosphorylate the downstream kinases.

Immunogen	A synthetic peptide of human CAMKK2 is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	Western Blot (Transfected lysate); ELISA
Size	1 ea
Buffer	In 1x PBS, pH 7.4
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	CAMKK2 calcium/calmodulin-dependent protein kinase kinase 2, beta [Homo sapiens]
Official Symbol	CAMKK2
Synonyms	CAMKK2; calcium/calmodulin-dependent protein kinase kinase 2, beta; calcium/calmodulin-dependent protein kinase kinase 2; CAMKK; CAMKKB; KIAA0787; MGC15254; caM-KK 2; caM-

KK beta; CAMKK beta protein; caM-kinase kinase 2; caM-kinase kinase beta;
calcium/calmodulin-dependent protein kinase beta;

Entrez Gene ID	10645
mRNA Refseq	NM_006549
Protein Refseq	NP_006540
UniProt ID	A0A024RBQ0
Chromosome Location	12q24.2
Pathway	Adipocytokine signaling pathway, organism-specific biosystem; Adipocytokine signaling pathway, conserved biosystem;
Function	ATP binding; calcium ion binding; calmodulin binding; calmodulin-dependent protein kinase activity; nucleotide binding; protein tyrosine kinase activity;
