



Anti-CCDC28A monoclonal antibody (DCABH-10849)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene is located in a region close to the locus of the pseudogene of chemokine (C-C motif) receptor-like 1 on chromosome 6. The specific function of this gene has not yet been determined.
Immunogen	A synthetic peptide of human CCDC28A is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	Western Blot (Transfected lysate); ELISA
Size	1 ea
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	CCDC28A coiled-coil domain containing 28A [Homo sapiens]
Official Symbol	CCDC28A
Synonyms	CCDC28A; coiled-coil domain containing 28A; C6orf80, chromosome 6 open reading frame 80; coiled-coil domain-containing protein 28A; CCRL1AP; DKFZp586D0623; chemokine C-C motif receptor-like 1 adjacent; C6orf80; MGC131913;
Entrez Gene ID	25901

Protein Refseq [NP_056254](#)

UniProt ID [Q8IWP9](#)

Chromosome Location 6q23.1-q24.1
