



Anti-CMTM8 monoclonal antibody (DCABH-11064)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene belongs to the chemokine-like factor gene superfamily, a novel family that is similar to the chemokine and the transmembrane 4 superfamilies. This gene is one of several chemokine-like factor genes located in a cluster on chromosome 3. This gene is widely expressed in many tissues, but the exact function of the encoded protein is unknown.
Immunogen	A synthetic peptide of human CMTM8 is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	Western Blot (Transfected lysate); ELISA
Size	1 ea
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	CMTM8 CKLF-like MARVEL transmembrane domain containing 8 [Homo sapiens]
Official Symbol	CMTM8
Synonyms	CMTM8; CKLF-like MARVEL transmembrane domain containing 8; chemokine like factor super family 8 , chemokine like factor superfamily 8 , CKLFSF8; CKLF-like MARVEL transmembrane domain-containing protein 8; chemokine-like factor superfamily 8; chemokine-like factor super family 8; chemokine-like factor superfamily member 8; CKLFSF8; CKLFSF8-V2;
Entrez Gene ID	152189

Protein Refseq	NP_849199
UniProt ID	Q8IZV2
Chromosome Location	3p22.2
Function	cytokine activity;