



Anti-FAM126A monoclonal antibody (DCABH-11488)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	The protein encoded by this gene may play a part in the beta-catenin/Lef signaling pathway. Expression of this gene is down-regulated by beta-catenin. Defects in this gene are a cause of hypomyelination with congenital cataract (HCC).
Immunogen	A synthetic peptide of human FAM126A is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	Western Blot (Transfected lysate); ELISA
Size	1 ea
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	FAM126A family with sequence similarity 126, member A [Homo sapiens]
Official Symbol	FAM126A
Synonyms	FAM126A; family with sequence similarity 126, member A; hyccin; down regulated by Ctnnb1; a; DRCTNNB1A; HCC; HYCC1; down regulated by Ctnnb1, a; down-regulated by CTNNB1 protein A; HLD5;
Entrez Gene ID	84668

Protein Refseq	NP_115970
UniProt ID	A0A024RA06
Chromosome Location	7p15.3
Function	signal transducer activity;