



# Anti-NAP1L3 monoclonal antibody (DCABH-12509)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

**Antigen Description** This gene encodes a member of the nucleosome assembly protein (NAP) family. It is linked closely to a region of genes responsible for several X-linked mental retardation syndromes. Alternatively spliced transcript variants have been found for this gene, but the biological validity of the variants has not been determined.

**Immunogen** A synthetic peptide of human NAP1L3 is used for rabbit immunization.

**Isotype** IgG

**Source/Host** Rabbit

**Species Reactivity** Human

**Purification** Protein A

**Conjugate** Unconjugated

**Applications** Western Blot (Transfected lysate); ELISA

**Size** 1 ea

**Buffer** In 1x PBS, pH 7.4

**Preservative** None

**Storage** Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

## GENE INFORMATION

**Gene Name** [NAP1L3 nucleosome assembly protein 1-like 3 \[ Homo sapiens \]](#)

**Official Symbol** NAP1L3

**Synonyms** NAP1L3; nucleosome assembly protein 1-like 3; MB20; MGC26312; NPL3;

**Entrez Gene ID** [4675](#)

**Protein Refseq** [NP\\_004529](#)

**UniProt ID** [Q8IYV1](#)

**Chromosome Location** Xq21.3-q22