



Anti-COMMD2 monoclonal antibody (DCABH-15039)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May down-regulate activation of NF-kappa-B.
Immunogen	A synthetic peptide of human COMMD2 is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	Western Blot (Transfected lysate); ELISA
Size	1 ea
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	COMMD2 COMM domain containing 2 [Homo sapiens]
Official Symbol	COMMD2
Synonyms	COMMD2; COMM domain containing 2; COMM domain-containing protein 2; HSPC042; MGC57611;
Entrez Gene ID	51122
Protein Refseq	NP_057178
UniProt ID	Q86X83
Chromosome Location	3q25.1

Function

protein binding;
