



Anti-PACRGL monoclonal antibody (DCABH-16313)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide of human PACRGL is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	Western Blot (Transfected lysate); ELISA
Size	1 ea
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	PACRGL PARK2 co-regulated-like [Homo sapiens]
Official Symbol	PACRGL
Synonyms	PACRGL; PARK2 co-regulated-like; C4orf28, chromosome 4 open reading frame 28; PACRG-like protein; MGC29898; C4orf28;
Entrez Gene ID	133015
Protein Refseq	NP_001124199
UniProt ID	Q8N7B6
Chromosome Location	4p15.31
Function	binding;