



Anti-PACSIN1 monoclonal antibody (DCABH-16314)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May play a role in vesicle formation and transport.
Immunogen	A synthetic peptide of human PACSIN1 is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	Western Blot (Transfected lysate); ELISA
Size	1 ea
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	PACSIN1 protein kinase C and casein kinase substrate in neurons 1 [Homo sapiens]
Official Symbol	PACSIN1
Synonyms	PACSIN1; protein kinase C and casein kinase substrate in neurons 1; protein kinase C and casein kinase substrate in neurons protein 1; SDPI; syndapin I; KIAA1379;
Entrez Gene ID	29993
Protein Refseq	NP_001186512
UniProt ID	Q5TZC3
Chromosome Location	6p21.3

Function

cytoskeletal protein binding; protein kinase activity;
