



Anti-TMTC2 monoclonal antibody (DCABH-17167)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	The function of this protein remains unknown.
Immunogen	A synthetic peptide of human TMTC2 is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	Western Blot (Transfected lysate); ELISA
Size	1 ea
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	TMTC2 transmembrane and tetratricopeptide repeat containing 2 [Homo sapiens]
Official Symbol	TMTC2
Synonyms	TMTC2; transmembrane and tetratricopeptide repeat containing 2; transmembrane and TPR repeat-containing protein 2; DKFZp762A217; IBDBP1 transcript 2; IBDBP1, transcript 1; IBDBP1;
Entrez Gene ID	160335
Protein Refseq	NP_689801
UniProt ID	Q8N394
Chromosome Location	12q21.31

Function

binding;
