



# Anti-TRIM74 monoclonal antibody (DCABH-17218)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	A synthetic peptide of human TRIM74 is used for rabbit immunization.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Protein A
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	Western Blot (Transfected lysate); ELISA
<b>Size</b>	1 ea
<b>Buffer</b>	In 1x PBS, pH 7.4
<b>Preservative</b>	None
<b>Storage</b>	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

## GENE INFORMATION

<b>Gene Name</b>	<a href="#">TRIM74 tripartite motif containing 74 [ Homo sapiens ]</a>
<b>Official Symbol</b>	TRIM74
<b>Synonyms</b>	TRIM74; tripartite motif containing 74; TRIM50C, tripartite motif containing 50C , tripartite motif containing 74; tripartite motif-containing protein 74; MGC45440; tripartite motif-containing 74; tripartite motif-containing 50C; tripartite motif-containing protein 50C; TRIM50C;
<b>Entrez Gene ID</b>	<a href="#">378108</a>
<b>Protein Refseq</b>	<a href="#">NP_942150</a>
<b>UniProt ID</b>	<a href="#">Q86UV6</a>
<b>Chromosome Location</b>	7q11.23
<b>Function</b>	metal ion binding; zinc ion binding;