



Anti-KRT18P55 monoclonal antibody (DCABH-200584)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	KRT18P55 (keratin 18 pseudogene 55) is a pseudogene, and is affiliated with the lncRNA class.
Immunogen	A synthetic peptide of human FLJ40504 is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	WB, ELISA
Size	1 mg
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	KRT18P55 keratin 18 pseudogene 55 [Homo sapiens (human)]
Official Symbol	KRT18P55
Synonyms	KRT18P55; keratin 18 pseudogene 55;
Entrez Gene ID	284085
Chromosome Location	17q11.2