



# Anti-SAP30BP monoclonal antibody (DCABH-201424)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Antigen Description</b>	SAP30BP (SAP30 binding protein) is a protein-coding gene. Diseases associated with SAP30BP include obesity, and schizophrenia. GO annotations related to this gene include protein binding. Induces cell death. May act as a transcriptional corepressor of a gene related to cell survival. May be involved in the regulation of beta-2-microglobulin genes.
<b>Immunogen</b>	A synthetic peptide of human SAP30BP is used for rabbit immunization.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Protein A
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, ELISA
<b>Size</b>	1 mg
<b>Buffer</b>	In 1x PBS, pH 7.4
<b>Storage</b>	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

## GENE INFORMATION

<b>Gene Name</b>	<a href="#">SAP30BP SAP30 binding protein [ Homo sapiens (human) ]</a>
<b>Official Symbol</b>	SAP30BP
<b>Synonyms</b>	SAP30BP; SAP30 binding protein; HTRG; HTRP; HCNGP; SAP30-binding protein; HSV-1 binding; transcriptional regulator protein HCNGP;
<b>Entrez Gene ID</b>	<a href="#">29115</a>
<b>mRNA Refseq</b>	<a href="#">NM_013260</a>

<b>Protein Refseq</b>	<a href="#">NP_037392</a>
<b>MIM</b>	<a href="#">610218</a>
<b>UniProt ID</b>	<a href="#">Q9UHR5</a>
<b>Chromosome Location</b>	17q25.1
<b>Function</b>	protein binding;