



Anti-TMEM54 monoclonal antibody (DCABH-201620)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	TMEM54 (transmembrane protein 54) is a protein-coding gene. Diseases associated with TMEM54 include hard palate cancer, and cervical adenocarcinoma. GO annotations related to this gene include molecular_function. An important paralog of this gene is KRTCAP3.
Immunogen	A synthetic peptide of human TMEM54 is used for rabbit immunization.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A
Conjugate	Unconjugated
Applications	WB, ELISA
Size	1 mg
Buffer	In 1x PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	TMEM54 transmembrane protein 54 [Homo sapiens (human)]
Official Symbol	TMEM54
Synonyms	TMEM54; transmembrane protein 54; BCLP; CAC1; CAC-1; protein CAC-1; beta-casein-like protein;
Entrez Gene ID	113452
Protein Refseq	NP_277039

UniProt ID	Q969K7
Chromosome Location	1p35-p34
Function	molecular_function;